
Specifications

Electrical Requirements

Model 525D		
Cat No. 11-690-525D	120 V, 50/60Hz	
Cat No. 11-690-526D	240 V, 50/60 Hz	
Model 537D		
Cat No. 11-690-537D	120 V, 50/60 Hz	
Cat No. 11-690-538D	240 V, 50/60 Hz	
Model 550D		
Cat No. 11-690-550D	120 V, 50/60 Hz	
Cat No. 11-690-551D	240 V, 50/60 Hz	

Power Requirements Chamber Volumes

Model 525D	260 W	2.5 cu ft
Model 537D	260 W	3.8 cu ft
Model 550D	320 W	5.0 cu ft

Chamber Dimensions (W x D x H)

Model 525D	18 x 18 x 13.5 in
Model 537D	18 x 18 x 20 in
Model 550D	18 x 18 x 26.5 in
