



|  |  |
|--|--|
| <b>Overall dimensions (L x W x H):</b> | 17 x 11 x 7.75"<br>(43.2 x 27.9 x 19.7cm)                      |
| <b>Tray dimensions (L x W):</b>        | 11 x 7.75"<br>(27.9 x 19.7cm)                                  |
| <b>Electrical (50/60 Hz):</b>          | 120 volts, 5 amps, 450 watts<br>230 volts, 5 amps, 450 watts   |
| <b>Fuses:</b>                          | 5mm x 20mm, 5 amp quick acting                                 |
| <b>Temperature range:</b>              | ambient +5°C to 65°C   |
| <b>Temperature uniformity:</b>         | ±0.5°C at 37°C   |
| <b>Speed range:</b>                    | 100 to 1200rpm   |
| <b>Speed accuracy:</b>                 | ±2% of set speed up to 999rpm<br>±5% 1000 to 1200rpm           |
| <b>Timer:</b>                          | 1 second to 9999 minutes<br>(increased in 1 second increments) |
| <b>Orbit:</b>                          | 0.125" (3mm)   |
| <b>Capacity:</b>                       | 4 microplates or 2 Micro-Tube Racks                            |
| <b>Controls:</b>                       | see page 4   |
| <b>Tray material:</b>                  | aluminum   |
| <b>Ship weight:</b>                    | 30lbs (13.6kg)   |