

7 Reference information

About this chapter

This chapter contains technical data, regulatory and other information.

7.1 Specifications

Parameter	Value
Ingression protection	IP21 (Frac-950) IP22 (Frac-920)
Supply voltage	100-240 V AC \pm 10%, 50 to 60 Hz
Power consumption	300 VA (Frac-950) 20 VA (Frac-920)
Dimensions (H x W x D)	480 x 380 x 550 mm (Frac-950) 250 x 320 x 400 mm (Frac-920)
Weight	16.5 kg (Frac-950) 4.5 kg (Frac-920)
Ambient temperature	4° to 40°C
Relative humidity tolerance	20% to 95% (non-condensing)
Atmospheric pressure	84 to 106 kPa (840 to 1060 mbar)
Acoustic noise level	73 dB A