



Overview

Catalog Number: 7754030

Specifications

- **Weight:** 275.0 lbs
- **Weight metric:** 125.0 kg
- **Dimensions:** 32.2" w x 28.0" d x 36.3" h
- **Dimensions metric:** 81.8 cm w x 71.1 cm d x 92.2 cm h
- **Electrical:** 230 volts, 50 Hz, 9.0 A
- **Region:** International
- **Collector Temperature:** -50Â° C, -58Â° F
- **Conformance:** CE
- **Ice Holding Capacity:** 12L
- **Plug Type:** Schuko
- **Style:** Console
- **Vacuum Pump Displacement Required:** 98 L/min or larger displacement

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative

