

HICLAVE HVE-50

◎Main specifications

Model No.	HVE-50	
Product name	"HICLAVE" Autoclave	
Outer Dimensions (Height with lid open)	540W×1040H×530Dmm (1370Hmm)	
Chamber size (Capacity)	300φ×710Dmm (50ℓ)	
Power source	AC110V/120V/220V/230V/240V, 50/60Hz (Please specify voltage when ordering)	
Power consumption	2.0KW	
Net weight (approx.)	57kg	
Sterilization temperature range	105 to 135℃	
Maximum allowable pressure (Gauge pressure)	0.26Mpa	
Temperature display range	5 to 137℃	
Warming temperature range	45 to 60℃	
Chamber material	Stainless steel (SUS304)	
Sterilization timer	1 to 250 minutes	
Pressure gauge	0 to 0.4 Mpa	
Selective mode	Sterilization of agar (with Warming)	Preparation ⇨ Heating ⇨ Sterilization ⇨ Exhaust ⇨ Warming (pulse) ⇨ Completion
	Sterilization of liquids:	Preparation ⇨ Heating ⇨ Sterilization ⇨ Exhaust ⇨ Completion (pulse)
	Sterilization of solids/ medical instruments	Preparation ⇨ Heating ⇨ Sterilization ⇨ Exhaust ⇨ Completion
Safety devices/ Warning alarms	<ul style="list-style-type: none"> *Double sensing inter lock mechanism *Over-pressure power cut-off *Over-temperature power cut-off *Sterilizing time counting timer *Heater malfunction checker *Temperature sensor wire breakage detector *Securely closed lid check function *Exhaust bottle ready-to-use check system *Lack-of-water prevention device Electric leakage breaker and over current detector Pressure safety valve <p>(Note: The * marked functions show ERROR indication on the display when it works)</p>	
Standard accessories	Stainless steel wire baskets, Heater cover, Drain hose, Exhaust bottle, Caster stoppers.	
Optional accessories	Stainless steel wire baskets, Stainless steel buckets.	

※ Chart recorder is an option.

*Specifications subject to change without notice due to continuous improvement.



Manufacturer

HIRAYAMA Manufacturing Corp.

2-6-5 TOYONO-CHO, KASUKABE-SHI,
SAITAMA 344-0014, JAPAN
TEL:+81-48-735-1242 FAX:+81-48-733-2384

Dealer