

Specifications

Dimensions (W x D x H):	6.75 x 12.75 x 14.25 inches (171 x 324 x 362 mm)
Weight:	27 lbs/12 Kg
Sealing Temperature:	OFF, 100°C to 190°C (increments of 1.0°)
Temperature Accuracy:	+/- 1.0° C
Temperature Uniformity:	+/- 1.0° C
Sealing Time:	0.5 to 10 sec. (increments of 0.1 sec)
Compatible Plate Materials:	PP (Polypropylene) PS (Polystyrene)
Compatible Plate Types:	Standard Assay Plates Deep-Well Storage Plates PCR Plates: skirted, semi-skirted, non-skirted
Maximum plate height:	45mm
Compatible Sealing Film Types:	Foil-polypropylene laminate Clear polyester-polypropylene laminate Clear polymer Thin clear polymer Foil-laminate Foil
Electrical Requirements:	AC100V ~ 120V, 50/60Hz, 400 Watts, 3.15A AC220V ~ 240V, 50/60Hz, 400 Watts, 2A