

SPECIFICATIONS

POWER REQUIREMENTS:

Model 4630:	120 VAC, 50/60 Hz, 0.5 Amps, 60 Watts
Model 4630JPN:	100 VAC, 50/60 Hz, 0.6 Amps, 60 Watts
Model 4630-1(CE):	240 VAC, 50/60 Hz, 0.2 Amps, 50 Watts

NOTE: ACCEPTABLE SUPPLY VOLTAGE FLUCTUATION RANGE FOR PROPER EQUIPMENT OPERATION (E.G., $\pm 10\%$ OF NORMAL SUPPLY).

SHAKER MOTION: A three-dimensional action that combines rocking and circular motions with a maximum range of movement of 10-degrees from the horizontal.

SHAKING SPEED: 30 revolutions-per-minute.

TIMER: 60-minute timer adjustable in 1-minute intervals or for continuous mixing, use the **CONSTANT** setting.

DIMENSIONS:

PLATFORM:	6 $\frac{3}{4}$ " W x 5"D (17 x 13 cm)
OVERALL:	6-7/8" W x 8-1/8"D x 7"H (17 x 21 x 18 cm)

SHIPPING WEIGHT: 25 lbs. (11 kg)

UNIT'S ENVIRONMENTAL OPERATING CONDITIONS:

POLLUTION DEGREE:	2
INSTALLATION CATEGORY:	II
ALTITUDE:	2000 Meters MSL (Mean Sea Level)
HUMIDITY:	80% maximum, non-condensing
ELECTRICAL SUPPLY:	120VAC or 240VAC
VOLTAGE TOLERANCE:	$\pm 10\%$ of normal rated line
TEMPERATURE:	15°C to 40°C
PRODUCT USAGE:	This product is intended for use indoors only