

# FreeZone Plus 12 Liter Cascade Console Freeze Dry System



View online: <https://www.labconco.com/product/freezone-plus-12-liter-console-freeze-dry-system-2/312>



## Overview

---

Catalog Number: 7960040

## Specifications

---

- **Weight:** 319.0 lbs
- **Weight metric:** 145.0 kg
- **Dimensions:** 32.2" w x 28.0" d x 36.3" h
- **Dimensions metric:** 81.8 cm w x 71.1 cm d x 92.2 cm h
- **Electrical:** 230 volts, 60 Hz, 9.0 A
- **Region:** U.S. and Canada
- **Collector Temperature:** -84Â° C, -119Â° F
- **Conformance:** CAN/CSA C22.2, UL
- **Ice Holding Capacity:** 12L
- **Plug Type:** North America, 230V
- **Style:** Console
- **Vacuum Pump Displacement Required:** 98 L/min or larger displacement

## For Product Assistance

---



[Andy Miller](#)  
Senior Sales Representative



[Vivian Goldmon](#)  
Customer Service Representative



Labconco Corporation  
8811 Prospect Avenue  
Kansas City, MO 64132-2696  
Phone: (816) 333-8811

