



Overview

Compliance: UL 61010-1, CAN/CSA C22.2 No. 61010.1, ETL listed.

Catalog Number: 4420430


Specifications

- **Weight:** 239.0 lbs
- **Weight metric:** 108.0 kg
- **Dimensions:** 24.2"w x 27.5"d x 34.2-36.2"h
- **Dimensions metric:** 61.5 cm w x 69.9 cm d x 86.9-91.9 cm
- **Electrical:** 115 V, 60 Hz, 16 amps. 20 amp circuit breaker is required.
- **Style:** Freestanding
- **Conformance:** CAN/CSA C22.2, ETL, UL
- **Viewing Window:** Yes

For Product Assistance



[KeiShaun Johnson](#)
Customer Coordinator



[Ryan Caton](#)
Sales Representative

