

Appendix C: Freeze Dryer Specifications

Catalog Number	Description	Voltage		Frequency	Phase	Amperage	
		Nominal	Operating Range			w/o Vacuum Pump	Max with Vacuum Pump
7753042	6 Liter Freeze Dryer	230	187-253	60	1	3	8
7753044	6 Liter Freeze Dryer	230	187-253	60	1	3	8
7753046	6 Liter Freeze Dryer	230	187-253	60	1	3	8
7753047	6 Liter Freeze Dryer	230	187-253	60	1	3	8
7753060	6 Liter Freeze Dryer	230	198-254	50	1	3	8
7753061	6 Liter Freeze Dryer	230	198-254	50	1	3	8
7753062	6 Liter Freeze Dryer	230	198-254	50	1	3	8
7753064	6 Liter Freeze Dryer	230	198-254	50	1	3	8
7753066	6 Liter Freeze Dryer	230	198-254	50	1	3	8
7753067	6 Liter Freeze Dryer	230	198-254	50	1	3	8
7753070	6 Liter Freeze Dryer	230	198-254	50	1	3	8
7753071	6 Liter Freeze Dryer	230	198-254	50	1	3	8
7753072	6 Liter Freeze Dryer	230	198-254	50	1	3	8
7753074	6 Liter Freeze Dryer	230	198-254	50	1	3	8
7753076	6 Liter Freeze Dryer	230	198-254	50	1	3	8
7753077	6 Liter Freeze Dryer	230	198-254	50	1	3	8
7753522	6 Liter Freeze Dryer	115	103-127	60	1	13	16
		100	90-110	50	1	13	16
7753524	6 Liter Freeze Dryer	115	103-127	60	1	13	16
		100	90-110	50	1	13	16
7753532	6 Liter Freeze Dryer	230	198-254	50	1	5.5	10
7753534	6 Liter Freeze Dryer	230	198-254	50	1	5.5	10
7753542	6 Liter Freeze Dryer	230	187-253	60	1	5.5	10
7753544	6 Liter Freeze Dryer	230	187-253	60	1	5.5	10
7753562	6 Liter Freeze Dryer	230	198-254	50	1	5.5	10
7753564	6 Liter Freeze Dryer	230	198-254	50	1	5.5	10
7753572	6 Liter Freeze Dryer	230	198-254	50	1	5.5	10
7753574	6 Liter Freeze Dryer	230	198-254	50	1	5.5	10
7934020	6 Liter Cascade Freeze Dryer	115	103-127	60	1	11	16
		100	90-110	50	1	11	16
7934021	6 Liter Cascade Freeze Dryer	115	103-127	60	1	11	16
		100	90-110	50	1	11	16
7934022	6 Liter Cascade Freeze Dryer	115	103-127	60	1	11	16
		100	90-110	50	1	11	16
7934024	6 Liter Cascade Freeze Dryer	115	103-127	60	1	11	16
		100	90-110	50	1	11	16
7934026	6 Liter Cascade Freeze Dryer	115	103-127	60	1	11	16
		100	90-110	50	1	11	16
7934027	6 Liter Cascade Freeze Dryer	115	103-127	60	1	11	16
		100	90-110	50	1	11	16

Appendix C: Freeze Dryer Specifications

Catalog Number	Description	Voltage			Amperage		
		Nominal	Operating Range	Frequency	Phase	w/o Vacuum Pump	Max with Vacuum Pump
7934030	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934031	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934032	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934034	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934036	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934037	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934040	6 Liter Cascade Freeze Dryer	230	187-253	60	1	5.5	10
7934041	6 Liter Cascade Freeze Dryer	230	187-253	60	1	5.5	10
7934042	6 Liter Cascade Freeze Dryer	230	187-253	60	1	5.5	10
7934044	6 Liter Cascade Freeze Dryer	230	187-253	60	1	5.5	10
7934046	6 Liter Cascade Freeze Dryer	230	187-253	60	1	5.5	10
7934047	6 Liter Cascade Freeze Dryer	230	187-253	60	1	5.5	10
7934060	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934061	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934062	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934064	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934066	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934067	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934070	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934071	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934072	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934074	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934076	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7934077	6 Liter Cascade Freeze Dryer	230	198-254	50	1	5.5	10
7754030	12 Liter Freeze Dryer	230	198-254	50	1	4.5	9
7754031	12 Liter Freeze Dryer	230	198-254	50	1	4.5	9
7754032	12 Liter Freeze Dryer	230	198-254	50	1	4.5	9
7754034	12 Liter Freeze Dryer	230	198-254	50	1	4.5	9
7754036	12 Liter Freeze Dryer	230	198-254	50	1	4.5	9
7754037	12 Liter Freeze Dryer	230	198-254	50	1	4.5	9
7754040	12 Liter Freeze Dryer	230	187-253	60	1	5	9
7754041	12 Liter Freeze Dryer	230	187-253	60	1	5	9
7754042	12 Liter Freeze Dryer	230	187-253	60	1	5	9
7754044	12 Liter Freeze Dryer	230	187-253	60	1	5	9
7754046	12 Liter Freeze Dryer	230	187-253	60	1	5	9
7754047	12 Liter Freeze Dryer	230	187-253	60	1	5	9
7754060	12 Liter Freeze Dryer	230	198-254	50	1	4.5	9
7754061	12 Liter Freeze Dryer	230	198-254	50	1	4.5	9
7754062	12 Liter Freeze Dryer	230	198-254	50	1	4.5	9
7754064	12 Liter Freeze Dryer	230	198-254	50	1	4.5	9
7754066	12 Liter Freeze Dryer	230	198-254	50	1	4.5	9
7754067	12 Liter Freeze Dryer	230	198-254	50	1	4.5	9