

FreeZone 6 Liter Console Freeze Dry System



View online: <https://www.labconco.com/product/freezone-6-liter-console-freeze-dry-system-with-purge-valve/152>



Overview

Includes purge valve with switch for isolating the vacuum pump from the freeze dry system.

Catalog Number: 7753022


Specifications

- **Weight:** 257.0 lbs
- **Weight metric:** 117.0 kg
- **Dimensions:** 32.2" w x 28.0" d x 36.3" h
- **Dimensions metric:** 81.8 cm w x 71.1 cm d x 92.2 cm h
- **Electrical:** 115 volts, 60 Hz, 16A
- **Region:** U.S. and Canada
- **Collector Temperature:** -50° C, -58° F
- **Conformance:** CAN/CSA C22.2, UL
- **Ice Holding Capacity:** 6L
- **Options Included:** Purge Valve
- **Plug Type:** North America, 115V
- **Style:** Console
- **Vacuum Pump Displacement Required:** 98 L/min or larger displacement

For Product Assistance



[KeiShaun Johnson](#)
Customer Coordinator



[Ryan Caton](#)
Sales Representative



Labconco Corporation
8811 Prospect Avenue
Kansas City, MO 64132-2696
Phone: (816) 333-8811

