



Overview

Manifolds are made of type 304 stainless steel and mount atop FreeZone Freeze Dryers.

Manifold with four neoprene valves with plastic knobs that accommodate both 1/2" and 3/4" adapters for connection of flask. Includes one latching pin.

Catalog Number: 7201700

Specifications

- **Weight:** 5.5 lbs
- **Weight metric:** 2.5 kg
- **Dimensions:** 8.6" w x 8.6" d x 11.2" h
- **Dimensions metric:** 21.8 x 21.8 x 28.4 cm
- **Applicable Products:** FreeZone 1 Liter, FreeZone 12 Liter, FreeZone 18 Liter, FreeZone 2.5 Liter, FreeZone 4.5 Liter, FreeZone 6 Liter, FreeZone 8 Liter

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative

