



## Overview

Manifolds are made of type 304 stainless steel and mount atop FreeZone Freeze Dryers.

Manifold with 14 neoprene valves with plastic knobs that accommodate both 1/2" and 3/4" adapters for connection of flask. Includes one latching pin. For use with FreeZone Console Freeze Dryers only.

**Catalog Number:** 7201900

## Specifications

- **Weight:** 21.0 lbs
- **Weight metric:** 9.6 kg
- **Dimensions:** 24.3" w x 10.6" d x 10.3" h
- **Dimensions metric:** 61.7 x 26.9 x 26.2 cm
- **Applicable Products:** FreeZone 12 Liter, FreeZone 18 Liter, FreeZone 6 Liter

## For Product Assistance



[Andy Miller](#)  
Senior Sales Representative



[Vivian Goldmon](#)  
Customer Service Representative

