



Overview

Accessories for use with the FreeZone Stoppering Tray Dryer, FreeZone Bulk Tray Dryer and FreeZone Small Tray Dryer are offered.

6-Port Tray Dryer Manifold and Isolation Valve. Six valves on manifold provide the flexibility to connect flasks and other freeze dry glassware to a Freeze Dryer with Stoppering or Bulk Tray Dryer (sold separately). Valves accommodate either 1/2" or 3/4" flask adapters. Separate valve isolates Tray Dryer from Freeze Dryer. This allows the user to pull samples under vacuum faster, determine endpoint, defrost ice in the middle of a cycle, and troubleshoot leaks easier.

Catalog Number: 7204300

Mounts to support stand included with the Tray Dryer. May be field installed on new or existing Freeze Dryers. Includes two couplings, two clamps and instruction sheet.

Specifications

- **Weight:** 45.0 lbs
- **Weight metric:** 20.4 kg
- **Dimensions:** 26.0" w x 32.0" d x 22.0" h
- **Dimensions metric:** 66.0 x 81.3 x 55.9 cm
- **Applicable Products:** Bulk Tray Dryers, Stoppering Tray Dryers

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative

