



## Overview

---

Provide connection of ampules to drying chambers and manifolds.

Attaches to a single valve port and accommodates up to 15 ampules for simultaneous freeze drying. Polypropylene base has 15 each 3/8" ID tubing ports, each with silicone tubing for attachment of ampules. For ampule neck sizes 10mm or larger. Includes 15 plugs, 2 each clamps and 1 straight stainless steel adapter. Ampules not included. Adapter attaches to port on freeze dryer.

**Catalog Number:** 7205101

## Specifications

---

- **Weight:** 2.5 lbs
- **Weight metric:** 1.3 kg
- **Dimensions:** 3.8" w x 3.8" d x 7.5" h
- **Dimensions metric:** 9.7 x 9.7 x 19.1 cm
- **Applicable Products:** Ampules

## For Product Assistance

---



[Andy Miller](#)  
Senior Sales Representative



[Vivian Goldmon](#)  
Customer Service Representative

