

Overview



The all-in-one CentriVac Complete Vacuum Concentrator houses every component you need to begin processing multiple samples; vacuum concentrator, -50°C cold trap, diaphragm pump, rotor and CentriZap strobe light. Just add your samples in the containers of your choice and start concentrating. Using centrifugal force, vacuum and controlled heat, the CentriVac Complete speeds concentration of samples used in a variety of applications found in molecular biology, proteomics, genomics, cell biology, microbiology, biochemistry, drug discovery and analytical chemistry.

Acid-resistant models are available for corrosive chemical use. Other models feature Heat Boost, which surrounds the chamber wall, to further speed evaporation of high boiling point solvents.

Catalog Number: 7315022

Specifications

- **Weight:** 145.0 lbs
- **Weight metric:** 66.0 kg
- **Dimensions:** 26.1" w x 22.7" d x 16.9" h
- **Dimensions metric:** 66.3 x 57.7 x 42.9 cm
- **Electrical:** 115 volts, 60 Hz, 12 amps, Domestic
- **Region:** U.S. and Canada
- **300-Watt Heat Boost:** Included
- **CentriVac Optional Accessories:** Rotors
- **CentriVac Required Accessories:** Sample Tubes
- **Conformance:** ETL
- **Lid Type:** Glass
- **Receptacle and Plug Type:** North America, 115 volts, 15 amps
- **Style:** Benchtop
- **Type/Application:** Aqueous

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldman](#)
Customer Service Representative

