

Acid-Resistant CentriVap Complete Vacuum Concentrator with Heat Boost



View online: <https://www.labconco.com/product/acid-resistant-centrivap-complete-vacuum-concentrator-with-heat-boost/3683>

Overview



The all-in-one CentriVap Complete Vacuum Concentrator houses every component you need to begin processing multiple samples; vacuum concentrator, -50°C cold trap, diaphragm pump, rotor and CentriZap strobe light. Just add your samples in the containers of your choice and start concentrating. Using centrifugal force, vacuum and controlled heat, the CentriVap Complete speeds concentration of samples used in a variety of applications found in molecular biology, proteomics, genomics, cell biology, microbiology, biochemistry, drug discovery and analytical chemistry.

Acid-resistant models are available for corrosive chemical use. Other models feature Heat Boost, which surrounds the chamber wall, to further speed evaporation of high boiling point solvents.

Catalog Number: 7315023

Specifications

- **Weight:** 145.0 lbs
- **Weight metric:** 66.0 kg
- **Dimensions:** 26.1" w x 22.7" d x 16.9" h
- **Dimensions metric:** 66.3 x 57.7 x 42.9 cm
- **Electrical:** 115 volts, 60 Hz, 12 amps, Domestic
- **Region:** U.S. and Canada
- **300-Watt Heat Boost:** Included
- **CentriVap Optional Accessories:** Rotors
- **CentriVap Required Accessories:** Sample Tubes
- **Conformance:** ETL
- **Lid Type:** Glass
- **Receptacle and Plug Type:** North America, 115 volts, 15 amps
- **Style:** Benchtop
- **Type/Application:** Acids, Aqueous

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative



Labconco Corporation
8811 Prospect Avenue
Kansas City, MO 64132-2696
Phone: (816) 333-8811

