



## Overview

---

Drying Chambers provide a means to attach flasks to your freeze dryer or contain samples in bulk.

12.0" diameter acrylic chamber for all FreeZone Freeze Dryers. Provides visibility for bulk freeze drying. Accommodates Tall Product Shelves (7317200). Includes aluminum lid with eight neoprene valves with plastic knobs for use with both 1/2" and 3/4" adapters, cap for vacuum port and pass through for electrical cords. 8" diameter opening to collector on bottom of chamber. Will fit 10.12" or 14.0" diameter condensers.

**Catalog Number:** 7318800

## Specifications

---

- **Weight:** 25.0 lbs
- **Weight metric:** 12.0 kg
- **Dimensions:** 17.4" diameter x 23.0" h
- **Dimensions metric:** 44.2 x 58.4 cm
- **Applicable Products:** FreeZone 1 Liter, FreeZone 12 Liter, FreeZone 18 Liter, FreeZone 2.5 Liter, FreeZone 4.5 Liter, FreeZone 6 Liter, FreeZone 8 Liter

## For Product Assistance

---



[Andy Miller](#)  
Senior Sales Representative



[Vivian Goldman](#)  
Customer Service Representative

