



Overview

CentriZap Strobe Light enables the user to shine a strobe light into the CentriVap's chamber to easily view samples while the rotor is spinning. Includes mounting holster and hardware for attachment to the right side of the CentriVap.

Catalog Number: 7464300

Specifications

- **Weight:** 2.0 lbs
- **Weight metric:** 0.9 kg
- **Applicable Products:** All CentriVap Benchtop and Mobile System Concentrators

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative

