

170 ml Borosilicate Glass Tube PTFE-Coated Aluminum Block



View online: <https://www.labconco.com/product/170-ml-borosilicate-glass-tube-ptfe-coated-aluminum-block/2380>



Overview

Blocks fit inside RapidVap Evaporation Systems. Glassware is required (not included).

Block tube capacity: 8. Sample tube size: 170 ml borosilicate glass tube. Sample volume: 125 ml. Actual sample volume depends on tube shape, height and vortex speed.

Catalog Number: 7486400

Specifications

- **Weight:** 8.2 lbs
- **Weight metric:** 3.7 kg
- **Applicable Products:** RapidVap N2 Evaporation Systems, RapidVap Vacuum Evaporation System
- **Block Tube Capacity:** 8

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative



Labconco Corporation
8811 Prospect Avenue
Kansas City, MO 64132-2696
Phone: (816) 333-8811

