

12mm OD Tube Size PTFE-Coated Aluminum Block



View online: <https://www.labconco.com/product/12mm-od-tube-size-ptfe-coated-aluminum-block/2374>

Overview



Blocks fit inside RapidVap Evaporation Systems. Glassware is required (not included).

Block tube capacity: 110. Sample tube size: 12 mm OD tubes, up to 6 ml. Maximum tube height: 150 mm. Sample volume: 4.5 ml. Actual sample volume depends on tube shape, height and vortex speed.

Catalog Number: 7491300

Specifications

- **Weight:** 8.0 lbs
- **Weight metric:** 3.6 kg
- **Applicable Products:** RapidVap Vacuum Evaporation System
- **Block Tube Capacity:** 110

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative



Labconco Corporation
8811 Prospect Avenue
Kansas City, MO 64132-2696
Phone: (816) 333-8811

