



## Overview

Product Shelves accommodate bulk samples, microtiter plates, vials or serum bottles. They are designed to fit inside Drying Chambers (sold separately).

Heated Product Shelves consist of three aluminum shelves with 1/4" rims, three stainless steel support rods, nine clips for shelf height adjustment, low voltage transformer with variable heat control to +40° C (+104° F), and three rubber feet. For use with Chambers 7522800, 7522900, 7443500, 7444000, 7442900 and 7867000.

**Catalog Number:** 7509200

## Specifications

- **Weight:** 13.0 lbs
- **Weight metric:** 5.9 kg
- **Dimensions:** 9.0" h x 6.5" diameter
- **Electrical:** 115 volts, 50/60 Hz, Domestic
- **Region:** U.S. and Canada
- **Applicable Products:** 12-Port Drying Chamber, 16-Port Drying Chamber, Large Clear Chamber, Large Clear Chamber with Valves, Small Clear Chamber, Tall Clear Chamber with Valves

## For Product Assistance



[Andy Miller](#)  
Senior Sales Representative



[Vivian Goldmon](#)  
Customer Service Representative

