



## Overview

Drying Chambers provide a means to attach flasks to your freeze dryer or contain samples in bulk.

Heated Drying Chamber is stainless steel with a clear acrylic cover, neoprene gaskets and electrical passthrough for transformer cords. It includes a three shelf product heater. The chamber is compatible with any FreeZone 4.5, 6, 12 or 18 Liter Freeze Dry System manufactured before August 2017 or any FreeZone Legacy Freeze Dryer manufactured after August 2017. If used with a FreeZone 4.5 Liter Freeze Dry System with permanent drying chamber, Attachment Port Lid Accessory 7762800 is required (not included). For operation on 115V, 50/60 Hz.

**Catalog Number:** 7521000

## Specifications

- **Weight:** 19.0 lbs
- **Weight metric:** 8.6 kg
- **Dimensions:** 10.5" h x 9" diameter
- **Electrical:** 115 volts, 50/60 Hz, Domestic
- **Region:** U.S. and Canada
- **Applicable Products:** FreeZone 12 Liter, FreeZone 18 Liter, FreeZone 4.5 Liter, FreeZone 6 Liter

## For Product Assistance



[Andy Miller](#)  
Senior Sales Representative



[Vivian Goldmon](#)  
Customer Service Representative

