

16-Port Drying Chamber



View online: <https://www.labconco.com/product/16-port-drying-chamber/3128>



Overview

Drying Chambers provide a means to attach flasks to your freeze dryer or contain samples in bulk.

The 16-Port Drying Chamber includes 3/4" thick acrylic lid with neoprene gasket, 16 neoprene valves that accommodate both 1/2" and 3/4" adapters for connection of flasks, and one latching pin. Made of type 304 stainless steel. It is compatible with FreeZone 1, 2.5, 4.5, 6, 8, 12 or 18 Liter Freeze Dryers.

Catalog Number: 7522900

Specifications

- **Weight:** 21.0 lbs
- **Weight metric:** 9.5 kg
- **Dimensions:** 13.0" diameter x 13.0" h
- **Applicable Products:** Freeze Dryers

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldman](#)
Customer Service Representative



Labconco Corporation
8811 Prospect Avenue
Kansas City, MO 64132-2696
Phone: (816) 333-8811

