

16-Port Drying Chamber with PTFE



View online: <https://www.labconco.com/product/16-port-drying-chamber-ptfe/6574>



Overview

Drying Chambers provide a means to attach flasks to your freeze dryer or contain samples in bulk.

The 16-Port Drying Chamber with PTFE is made of type 304 stainless steel with a corrosion-resistant PTFE-coated interior. Includes 16 neoprene valves with plastic knobs that accommodate both 1/2" and 3/4" adapters for connection of flasks, a 3/4" thick acrylic lid with neoprene gasket, and one latching pin. It is compatible with FreeZone 1, 2.5, 4.5, 6, 8, 12 or 18 Liter Freeze Dryers.

Catalog Number: 7522901

Specifications

- **Weight:** 21.0 lbs
- **Weight metric:** 10.0 kg
- **Dimensions:** 13.0" diameter x 13.0" h
- **Dimensions metric:** 33.0 x 33.0 cm
- **Applicable Products:** Freeze Dryers

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative



Labconco Corporation
8811 Prospect Avenue
Kansas City, MO 64132-2696
Phone: (816) 333-8811

