



Overview

Drying Chambers provide a means to attach flasks to your freeze dryer or contain samples in bulk.

The Mini Stoppering Chamber attaches to a port on a Labconco Drying Chamber or Manifold (sold separately) to stopper small vials (up to 30ml each). Includes 17" of clear PVC tubing.

Can stopper up to 48 serum vials (2ml) or 8 serum vials (30ml) while maintaining vacuum. Spacers, tubing, and adapter for connection to the port are included. Vials and stoppers sold separately.

Catalog Number: 7569000

Specifications

- **Weight:** 15.0 lbs
- **Weight metric:** 6.8 kg
- **Dimensions:** 12" w x 12" d x 10" h
- **Dimensions metric:** 30.5cm x 30.5cm x 25.4cm
- **Applicable Products:** Drying Chambers, Manifolds

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative

