


Overview



Accessories for use with the FreeZone Stoppering Tray Dryer, FreeZone Bulk Tray Dryer and FreeZone Small Tray Dryer are offered.

6-Port Tray Dryer Manifold. Six valves provide the flexibility to connect flasks and other freeze dry glassware to the Freeze Dryer with Stoppering or Bulk Tray Dryer. Valves accommodate either 1/2" or 3/4" flask adapters. Mounts to support stand included with the Tray Dryer. May be field installed before or after the Tray Dryer is installed on the Freeze Dryer. Included with Bulk Tray Dryer 7806021 and 7806031.

Catalog Number: 7726500

Specifications

- **Weight:** 15.0 lbs
- **Weight metric:** 6.8 kg
- **Dimensions:** 30" w x 16.5"d x 7" h
- **Dimensions metric:** 76.2 x 41.9 cm
- **Applicable Products:** Bulk Tray Dryers, Stoppering Tray Dryers

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative

