



Overview

For use with any FreeZone Freeze Dryer manufactured after July 2017. Chambers allow for drying of bulk samples, or drying of serum bottles, vials or microplates. The 10" diameter shelves are wireless and removable so that they may be placed in a freezer for pre-freezing. Each shelf is height adjustable to accommodate a broad variety of sample sizes. Chamber connects to a FreeZone Freeze Dryer for control and recording of shelf temperature from the Lyo-Works display. Temperature may be set from -20° C (-4° F) to +60° C (+140° F).

Each Chamber includes a 12" diameter clear acrylic cylinder, vacuum release valve, cap for vacuum port and 8" diameter opening to collector on bottom of chamber. Will fit 10.12" or 14.0" diameter condensers. Short Chambers include three shelves and three temperature probes. Tall Chambers include five shelves and five temperature probes.

Catalog Number: 780801000

Specifications

- **Weight:** 80.0 lbs
- **Weight metric:** 36.0 kg
- **Dimensions:** 13.5" w x 15.2" d x 17.0" h
- **Dimensions metric:** 34.3 x 38.6 x 43.2 cm
- **Electrical:** 115V, 50/60 Hz, 2A
- **Region:** U.S. and Canada
- **Applicable Products:** FreeZone 12 Liter, FreeZone 18 Liter, FreeZone 2.5 Liter, FreeZone 4.5 Liter, FreeZone 6 Liter, FreeZone 8 Liter
- **Conformance:** CAN/CSA C22.2, UL
- **Plug Type:** North America, 115V

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldman](#)
Customer Service Representative

