



Overview

FreeZone Benchtop Shell Freezers provide a means to freeze a thin coating of sample on the inside walls of flasks. Since shell frozen samples have a greater surface area than slab frozen samples, freeze drying is enhanced.

Catalog Number: 7949020

Specifications

- **Weight:** 95.0 lbs
- **Weight metric:** 43.1 kg
- **Dimensions:** 14.5" w x 22.1" d x 19.7" h
- **Dimensions metric:** 36.8 x 56.1 x 50.0 cm
- **Electrical:** 115 volts, 60 Hz, 6 A
- **Region:** U.S. and Canada

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative

