



## Overview

Vacuum Controller provides manual control and monitoring of vacuum levels. The vacuum controller features an analog gauge shows vacuum from 0-30" Hg. At the end of concentration, the control valve may be opened to break the vacuum. For use with Diaphragm Vacuum Pumps only.

**Catalog Number:** 7991000

## Specifications

- **Weight:** 3.0 lbs
- **Weight metric:** 1.4 kg
- **Applicable Products:** Diaphragm Vacuum Pumps

## For Product Assistance



[Andy Miller](#)  
Senior Sales Representative



[Vivian Goldman](#)  
Customer Service Representative

