

SPECIFICATIONS

POWER REQUIREMENTS:

| | |
|-------------------------|--|
| Model 4625: | 120 VAC, 60 Hz, 120 Watts, 1Amp |
| 4625JPN: | 100 VAC, 50/50 Hz, 120 Watts, 1.2 Amps |
| Model 4625-1, 4625-1CE: | 240 VAC, 50 Hz, 120 Watts, 0.5 Amp |
| Model 4625-2: | 240 VAC, 60 Hz, 120 Watts, 0.5 Amp |

MOTOR: Continuous-duty rated, high-efficiency DC motor features high torque, full ball-bearing construction, and ceramic permanent magnets. Brushes are user replaceable.

SHAKING SPEED: Speed range is 40 to 1100 rpm.

TIMER: Range: Up to 5 minutes in 1 minute increments or continuous run operation.

SHAKER MOTION: 0.120-inch orbit

DIMENSIONS: 11" W x 11" D x 7" H (28 x 28 x 18 cm)

SHIPPING WEIGHT: 25 lbs. (11 kg.)

UNIT'S ENVIRONMENTAL OPERATING CONDITIONS:

| | |
|------------------------|--|
| POLLUTION DEGREE: | 2 |
| INSTALLATION CATEGORY: | II |
| ALTITUDE: | 2000 Meters MSL (Mean Sea Level) |
| HUMIDITY: | 80% maximum, non-condensing |
| ELECTRICAL SUPPLY: | 120VAC or 240VAC |
| VOLTAGE TOLERANCE: | ±10% of normal rated line |
| TEMPERATURE: | 15°C to 40°C |
| PRODUCT USAGE: | This product is intended for use indoors only |