

## Appendix D

### Product specifications

Feature	Specifications
Input Power (switchable)	115 VAC, 50/60 Hz 230 VAC, 50/60 Hz
Fuses	One 4A/250V, one extra fuse is provided
Output power in watts	300 watts
Output voltage range	5~250 V
Output current range	10 mA~3 A
Timer	~99.99 per step
Terminal pairs	4 (4 positive voltage and 4 negative voltage)
Operating Modes <ul style="list-style-type: none"> <li>• Constant Voltage</li> <li>• Constant Current</li> <li>• Constant Power</li> </ul>	1V step 10 mA step 1 watt step
Crossover	Automatic
Display type	Backlit LCD Graphic type
Display size	53.64 x 15.64 mm (W x H)
Pause function	Yes
Safety features	No Load Detection Load Change Detection Overload Detection Ground Leak Detection Auto Restart
Programmable	Yes <ul style="list-style-type: none"> <li>• Store file no.: 20</li> <li>• Program: up to 10 steps</li> </ul>
Stackable	Yes
Housing material	Flame retardant ABS
Housing size	200 x 290 x 70 mm (W x D x H)
Operating temperature	0°C–40°C
Environmental condition	≤95% RH (non-condensing), 75 KPa–106 Kpa, altitude not to exceed 2000 meters
Weight	2.6 kg
Certifications	CE, cETLus, c-Tick
Warranty	1 year