

Cat. No.	Capacity	Interior D x W x H	Exterior L x W x H	Temp. Resolution/Uniformity	Temp. Range	No. of Shelves	Finish	Electrical	Shipping Weight
VO914A	18.6L (0.65 cu. ft.)	36 × 23 × 23cm (14 × 9 × 9in.)	58 × 46 × 66cm (23 × 18 × 26in.)	±1°C/ ±3.5°C	6° above ambient to 260°C	3	Painted Steel	120V 50/60Hz 750w	68kg (150 lb.)
VO914C	18.6L (0.65 cu. ft.)	36 × 23 × 23cm (14 × 9 × 9in.)	58 × 46 × 66cm (23 × 18 × 26in.)	±1°C/ ±3.5°C	6° above ambient to 260°C	3	Painted Steel	208/240V 50/60Hz 750w	68kg (150 lb.)
VO914SA	18.6L (0.65 cu. ft.)	36 × 23 × 23cm (14 × 9 × 9in.)	58 × 46 × 66cm (23 × 18 × 26in.)	±1°C/ ±3.5°C	6° above ambient to 260°C	3	Stainless Steel	120V 50/60Hz 750w	68kg (150 lb.)
VO1218A	42.5L (1.5 cu. ft.)	46 × 31 × 31cm (18 × 12 × 12in.)	66 × 53 × 74cm (26 × 21 × 29in.)	±1°C/ ±3.5°C	6° above ambient to 260°C	3	Painted Steel	120V 50/60Hz, 1250w	113kg (250 lb.)
VO1218C	42.5L (1.5 cu. ft.)	46 × 31 × 31cm (18 × 12 × 12in.)	66 × 53 × 74cm (26 × 21 × 29in.)	±1°C/ ±3.5°C	6° above ambient to 260°C	3	Painted Steel	208/240V 50/60Hz 1250w	113kg (250 lb.)
VO1218SA	42.5L (1.5 cu. ft.)	46 × 31 × 31cm (18 × 12 × 12in.)	66 × 53 × 74cm (26 × 21 × 29in.)	±1°C/ ±3.5°C	6° above ambient to 260°C	3	Stainless Steel	120V 50/60Hz, 1250w	113kg (250 lb.)
VO1824A	127.4L (4.5 cu. ft.)	61 × 46 × 46cm (24 × 18 × 18in.)	81 × 69 × 89cm (32 × 27 × 35in.)	±1°C/ ±3.5°C	6° above ambient to 260°C	2	Painted Steel	120V 50/60Hz 1500w	147kg (325 lb.)
VO1824C	127.4L (4.5 cu. ft.)	61 × 46 × 46cm (24 × 18 × 18in.)	81 × 69 × 89cm (32 × 27 × 35in.)	±1°C/ ±3.5°C	6° above ambient to 260°C	2	Painted Steel	208/240V 50/60Hz 1500w	147kg (325 lb.)
VO1824SA	18.6L (0.65 cu. ft.)	61 × 46 × 46cm (24 × 18 × 18in.)	81 × 69 × 89cm (32 × 27 × 35in.)	±1°C/ ±3.5°C	6° above ambient to 260°C	2	Stainless Steel	120V 50/60Hz, 1500w	147kg (325 lb.)
VO1824HPC	127.4L (4.5 cu. ft.)	61 × 46 × 46cm (24 × 18 × 18in.)	81 × 69 × 89cm (32 × 27 × 35in.)	±1°C/ ±3.5°C	6° above ambient to 260°C	1	Painted Steel	208/240V 50/60Hz 3000w	147kg (325 lb.)

Thermo Scientific Accessories for Lindberg/Blue M Vacuum Ovens

Catalog No	Description	For use with
6718	White Floor Stand	All Lindberg/Blue M vacuum ovens
6723	Stainless-steel Floor Stand	All Lindberg/Blue M vacuum ovens
118961	Connection Kit – includes NW fittings, centering rings, connectors, vacuum clamps and hose to connect vacuum pump to oven.	All Lindberg/Blue M vacuum ovens
305994H01	Silicone replacement gasket (red)	VO914 Series
305992H01	Silicone replacement gasket (red)	VO1218 Series
305993H01	Silicone replacement gasket (red)	VO1824 Series
118974	Buna-N silicone-free gasket, for temperatures up to 150°C only	VO914 Series
118977	Buna-N silicone-free gasket, for temperatures up to 150°C only	VO1218 Series
34637H01	Buna-N silicone-free gasket, for temperatures up to 150°C only	VO1824 Series
6826	Vacuum Pump	For VO914 Series
6827	Vacuum Pump	For VO1218 and 1824 Series
6831	Vacuum Pump Oil	All Vacuum pumps
6832	Mist Filter Element	All Vacuum pumps
6834	Activated Alumina	All Lindberg/Blue M vacuum ovens