

## 11.Specifications

Items	Specifications
Voltage [VAC]	100-120/200-240
Frequency [Hz]	50/60
Power [W]	1010
Plate material	Glass ceramic
Dimensions of workplate (mm)	184×184
*Heating power [W]	1000
*Temperature range [ °C]	RT-550, increment: 5
*Temperature display [ °C]	LED
*Temperature display accuracy [ °C]	± 1
*Control accuracy of heating temperature [ °C]	± 10
*The safety temperature of the hotplate [ °C]	580

*Temperature sensor in medium	PT1000
*Control accuracy of heating temperature with temperature sensor [ °C]	± 0.5
*Residual heat warning	50°C
Dimensions (mm)	215×360×112
Weight [kg]	3.8
Permitted ambient temperature [ °C]	5-40
Permitted relative humidity	80%
Protection class acc. to DIN 60529	IP21

Table 3