

## Section 3 Specifications

### Electrical Requirements

SHKA2506 120 Volt, 50/60 Hertz, 40 Watts, 0.45 Amp

SHKA2506-1 240 Volt, 50/60 Hertz, 40 Watts, 0.45 Amp

**Gear Motor** 1/18th HP, DC permanent magnet motor

**Shaking Movement** Reciprocating motion with adjustable stroke range in 1/4" increments from 1/2" to 1-1/2"

**Shaking Speed** 10-250 motions per minute

**Shaking Interval** Continuous or timed intervals from 1 minute to 60 minutes

**Dimensions** 13.5" W X 17.4" D X 7.2" H (34.3 X 44.2 x 18.3 cm)

**Shipping Weight** 27 lbs. (12.3 kg)

### Environmental Operating Conditions

Pollution Degree\*\* .....2

Installation Category\*\* .....II

Altitude .....2000 meters MSL (Mean Sea Level)

Relative Humidity . . .20% to 80% maximum, non-condensing

Electrical Supply .....120 VAC or 240 VAC

Voltage Tolerance .....+10% of normal rated line

Temperature .....0°C to 40°C

Product Usage . . . . This product is intended for use indoors only

*\*\*Refer to IEC 664-1*