

PAJ Gold Series Jewelry Balances

Models	PAJ812N	PAJ2102N	PAJ3102N	PAJ4102N	PA4101N
Capacity	810g	2100g	3100g	4100g	4100g
Readability	0.01g				0.1g
Certified Readability	0.1g				
Repeatability (Std Dev)	10mg			100mg	
Linearity	± 30mg			± 100mg	
Application Modes	Weighing, parts counting (with APW optimization)				
Weighing Units	g, lb, ct, oz, dwt, tical, tola, mommes, bhat, grain, mesghal, Newton, ozt		g, kg, lb, ct, oz, dwt, tical, tola, mommes, bhat, grain, mesghal, Newton, ozt		
Platform Size	6.6 x 7.1 in / 16.8 x 18 cm				
Calibration Weight-Span	400 / 800g	1000g / 2000g	2 / 3kg	2 / 4kg	2 / 4kg
Calibration Weight-Linearity	400, 800g	1000, 2000g	2, 3kg	2, 4kg	2, 4kg
Stabilization Time	Within 3 seconds				
Interface	RS232, bi-directional				
Power	AC adapter (included)				
Operating Conditions	14° F / -10° C to 104° F / 40° C at 10% to 80% relative humidity, non-condensing				
Storage Temperature	-40° F / -40° C to 158° F / 70° C at 10% to 80% relative humidity, non-condensing				
Net Weight	7.3 lb / 3.3 kg				
Shipping Weight	11.8 lb / 5.4 kg				
Balance Dimensions (W x H x D)	7.7 x 3.6 x 12.6 in / 19.6 x 9.2 x 32 cm				
Carton Dimensions (L x W x H)	19.5 x 15.5 x 12.8 in / 49.5 x 39.5 x 32.7 cm				

Other Standard Features and Equipment

AC adapter, user selectable span calibration points, software reset menu, stability indicator, Auto Tare, user selectable communications settings, user selectable printing options, stainless steel platform, weigh-below hook, integral security bracke

Compliance

Legal for Trade: NTEP Class II, CoC No. 09-043

Safety: FCC (Part 15 Class A), CE (EN 61326 Class A, Criteria A) C Tick N13123 (AS/NZS4251.1, AS/NZS4252.1)

Accessories

Dot Matrix Printer	80252042
Cable kit to connect Dot Matrix Printer to PAJ	80252571
Paper for Dot Matrix Printer	80251932
Ink Ribbon for Dot Matrix Printer	80251933
Thermal Printer	80251992
Paper for Thermal Printer	80251931
Cable kit to connect Thermal Printer to PAJ	80252581
OHAUS Collect Software	80500746
TAL Software WinWedge	SW12W
RS232 Cable, PC 9 Pin to PAJ	80500525
RS232 Cable, PC 25 Pin to PAJ	80500524
Auxiliary Display	30472064
Security Device for PAJ	76288-01

Industry Leading Quality and Support

All OHAUS PAJ balances are manufactured under an ISO 9001:2015 Registered Quality Management System. Our rugged construction and stringent quality control have been hallmarks of all OHAUS products for over a century.

OHAUS CORPORATION

7 Campus Drive
Suite 310
Parsippany, NJ 07054 USA

Tel: 800.672.7722
973.377.9000
Fax: 973.944.7177

www.ohaus.com

The management system governing the manufacture of this product is ISO 9001:2015 certified.

