

SK-R330-PRO/SK-D3307-PRO/SK-D3309-PRO

User Manual

SK-R330-Pro LCD digital Rocking Shaker

SK-D3307-Pro LCD digital 3D Shaker

SK-D3309-Pro LCD digital 3D Shaker



Please read the User Manual carefully before use, and follow all operating and safety instructions!

www.scilogex.com

SCIOGEX

9. Specifications

Items	Specifications	
	Rocking shaker	3D shaker
Voltage [VAC]	100-240	100-240
Frequency [Hz]	50/60	50/60
Power [W]	40	50
Tile Angel [°]	±9	±7, ±9
Max. Shaking Weight [kg]	7	5
Motor Type	DC motor	DC motor
Speed range [rpm]	10-70	10-70
Speed display	LCD	LCD
Timer range [min]	1-1199	1-1199
Timer display	LCD	LCD
Run type	Timer/Continuous operation	Timer/Continuous operation

Items	Specifications	
	Rocking shaker	3D shaker
Overall Dimensions [mm]	360×420×155	430×360×106
Weight [kg]	12	9
Permissible ambient temperature [°C]	5-40	5-40
Permissible relative humidity	80%	80%
Protection class acc. to DIN EN60529	IP21	IP21

Table 4