




ROTOR SPECIFICATIONS

Type	Swinging Bucket
Material	Aluminum
Capacity	4 x 750 mL
Maximum Speed	4,700 rpm
Tube Dimensions Ø x L	98 x 135 mm
Angle	90°
K Factor	8,138
Net Weight	6.0 kg/13.2 lbs
Warranty	10 years
Certifications¹	

	RCF	RADIUS
Maximum	4,575 xg (Buckets) 3,706 xg ² (Microplate Carriers)	185.4 mm (Buckets) 150.2 mm (Microplate Carriers)
Minimum	3,411 xg (Buckets)	91.1 mm (Buckets)

ROTOR AND ACCESSORIES ORDERING INFORMATION

SH-3000 Rotor Package includes:

Description	Cat. No.
SH-3000 Rotor	11796
750 mL Buckets (qty. 2)	11788
Rotor Extractor	01031
Instruction Manual	11782

Accessories:

Description	Cat. No.
Microplate Carriers (qty. 2)	11779
Sealing Covers for 750 mL Buckets ³	11776
Sealing Cover O-rings (qty. 6)	12995

LABWARE ORDERING INFORMATION

Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF (xg)	Dims ØxL (mm)	REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
									Type	Cat. No.	Qty. /Pkg.	Description
750.0	750.0	PP Bottle	75003699	12	4	4,700	4,575	98 x 135	Closure	Included	12	PP Sealing
500.0	—	Conical Tissue Culture Bottle	—	—	4	—	—	96 x 146	Adapter	Available from bottle manufacturer	1	1 place/adapter
250.0	250.0	PP Oak Ridge Bottle	3141-0250	36	4	4,700	4,575	61 x 133	Closure Adapter	Included 01092	36 1	PP Sealing 1 place/adapter
250.0	250.0	PC Oak Ridge Bottle	3140-0250	36	4	4,700	4,575	62 x 135	Closure Adapter	Included 01092	36 1	PP Sealing 1 place/adapter

¹  Biocontainment certification by the Public Health Laboratory Service, Centre for Applied Microbiological Research, Porton Down, UK, shown to be suitable for use with materials up to ACDP Category 3 as categorized by the Advisory Committee on Dangerous Pathogens.

² Maximum g-force rating of microplate varies depending upon manufacturer. Verify correct maximum speed/g-force before use.

³ The rotor/cover assembly has been tested for microbiological containment by PHLS-CAMR-Porton Down. Suitable for use with materials up to ACDP category 3. Seal integrity may be affected by improper use or maintenance.

LABWARE ORDERING INFORMATION CONTINUED

Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF (xg)	Dims ØxL (mm)	REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
									Type	Cat. No.	Qty. /Pkg.	Description
250.0	280.0	PP Dry-Spin Bottle	03069	12	4	4,700	4,575	61 x 146	Closure	Included	4	PP Sealing
									Closure	06841	12	Replacement Caps
									Closure	03545	24	Replacement O-Rings
250.0	270.0	PC Dry-Spin Bottle	03098	12	4	4,700	4,575	61 x 146	Adapter	01092	1	1 place/adapter
									Closure	Included	4	PP Sealing
									Adapter	01092	1	1 place/adapter
									Closure ⁴	06841	12	Replacement Caps
200.0	200.0	Nunc Tissue Culture Conical Bottle	03094	4	4	4,700	4,575	60 x 135	Closure	03545	24	Replacement O-Rings
									Adapter	01092	1	1 place/adapter
									Adapter	03095	1	1 place/adapter
									Sealing	Included	4	PP Non Sealing
175.0	175.0	PP Conical Culture Bottle	3143-0175	36	4	4,700	4,575	62 x 144	Adapter	—	1	1 place/adapter
									Adapter	03095	1	1 place/adapter
175.0	175.0	PC Conical Culture Bottle	3144-0175	36	4	4,700	4,575	62 x 144	Closure	DS3126-0175	1	1 place/adapter
									Adapter	DS3126-0175	1	1 place/adapter
175.0	175.0	PS Sterile Conical Culture Bottle	3145-0175	12	4	4,700	4,575	62 x 144	Closure	Included	4	PP Sealing
									Adapter	01092	1	1 place/adapter
									Adapter	DS3126-0175	1	1 place/adapter
150.0	144.0	PC Thick-Walled Tube	03143	25	4	4,700	4,575	45 x 132	Adapter	01092	1	1 place/adapter
									Adapter	00458	1	1 place/adapter
100.0	100.0	GL Pyrex® Tube	03230	1	4	3,000	1,864	44 x 137	Closure	Included	—	PP Non-Sealing
									Adapter ⁵	00371	1	1 place/adapter
									Adapter ⁵	01092	1	1 place/adapter
50.0	—	Conical Tissue	—	—	12	—	—	30 x 115	Adapter ⁶	00436	1	5 places/adapter
									Adapter ⁶	00445	1	3 places/adapter
50.0	37.0	PP Flanged Tube	03147	25	20	4,000	3,313	29 x 102	Closure	03268	25	PP Snap-on
									Adapter	00894	1	5 places/adapter Non-Sealing
50.0	36.0	PC Flanged Tube	03146	25	20	4,000	3,313	29 x 102	Closure	03268	25	PP Snap-on
									Adapter	00894	1	5 places/adapter
50.0	46.0	PA Thin-Walled Tube	03139	25	20	4,700	4,575	28 x 104	Adapter	00894	1	5 places/adapter
									Closure	Included	10	FEP Sealing
50.0	46.0	FEP Oak Ridge Tube	3114-0050	10	20	4,700	4,575	29 x 107	Adapter	00894	1	5 places/adapter
									Closure	Included	100	PP Non-Sealing
50.0	50.0	PP Oak Ridge Tube	3119-0050	100	20	4,700	4,575	29 x 107	Adapter	00894	1	5 places/adapter
									Closure	Included	100	PP Non-Sealing
50.0	50.0	PC Oak Ridge Tube	3118-0050	100	20	4,700	4,575	61 x 146	Adapter	00894	1	5 places/adapter
									Closure	Included	50	PP Sealing
50.0	43.0	PC Nalgene Oak Ridge Tube	3138-0050	50	20	4,700	4,575	29 x 107	Closure	Included	50	PP Sealing
									Closure	Included	50	PP Sealing
50.0	50.0	SS Flanged Tube	00517	1	20	4,700	4,575	28 x 101	Closure	00518	1	Stainless Steel Sealing
									Adapter	00894	1	5 places/adapter
									Tool	01014	1	Wrench
16.0	13.0	PP Flanged Tube	03244	50	48	4,700	—	18 x 100	Closure	03299	50	PP Snap-On
									Adapter	00892	1	12 places/adapter
16.0	14.0	PC Flanged Tube	03243	50	48	4,700	—	18 x 100	Closure	03299	50	PP Snap-On
									Adapter	00892	1	12 places/adapter
15.0	—	Conical Tissue Culture Tube	—	—	32	—	—	17 x 120	Adapter	00447	1	8 places/adapter
									Adapter	00385	1	18 places/adapter
14.0	11.0	PP Flanged Tube	03247	50	20	4,700	4,575	18 x 75	Closure	03269	50	PP Snap-On
									Adapter	00894	1	5 places/adapter
									Adapter	00382 ⁷	1	—

⁴ Replacement part. Not required for use.⁵ Both adapters must be used together.⁶ Both adapters will accommodate skirted and standard conical tissue culture tubes. With sealed covers, use adapter cat. no. 00445.

LABWARE ORDERING INFORMATION CONTINUED

Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF (xg)	Dims ØxL (mm)	REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
									Type	Cat. No.	Qty. /Pkg.	Description
14.0	11.0	PC Flanged Tube	03246	50	48	4,700	4,575	18 x 75	Closure	03269	50	PP Snap-On
									Adapter	00892	1	12 places/adapter
12.0	12.0	GL Pyrex Tube	03109	50	20	4,700	4,575	18x102	Adapter	00894	1	5 places/adapter
										00382	5	—
12.0	10.0	PP Flanged Tube	03116	50	72	4,700	4,575	16 x 100	Closure	03266	50	PP Nonsealing
									Adapter	00385	1	18 places/adapter
12.0	11.0	PC Flanged Tube	03115	50	72	4,700	4,575	16 x 100	Closure	03266	50	PP Nonsealing
									Adapter	00385	1	18 places/adapter
10.0	—	Blood Collection Tube	—	—	72	—	—	16 x 100	Adapter	00385	1	18 places/adapter
10.0	9.1	PP Oak Ridge Tube	03929	50	20	4,700	4,575	16 x 80	Closure	03279	25	PP Sealing
									Adapter	00894	1	5 places/adapter
										00402	5	—
10.0	9.1	PP Oak Ridge Tube	03020	25	20	4,700	4,575	16 x 80	Closure	03279	25	PP Sealing
									Adapter	00894	1	5 places/adapter
										00402	5	—
10.0	—	Blood Collection Tube	—	—	88	—	—	13 x 100	Adapter	00384	1	22 places/adapter
7.0	6.0	PP Flanged Tube	03121	50	88	4,700	4,575	13 x 100	Closure	03265	50	PP Snap-On
									Adapter	00384	1	22 places/adapter
7.0	6.0	PC Flanged Tube	03120	50	88	4,700	4,575	13 x 100	Closure	03265	50	PP Snap-On
									Adapter	00384	1	22 places/adapter
5.0	—	GL Standard Tube	—	—	132	—	—	12 x 75	Adapter	00857	1	33 places/adapter
5.0	—	Blood Collection Tube	—	—	132	—	—	13 x 75	Adapter	00857	1	33 places/adapter
5.0	5.0	GL Pyrex Tube	03102	50	132	4,700	4,575	12 x 75	Adapter	00857	1	33 places/adapter
4.0	—	GL Standard Glass Tube	—	—	132	—	—	10 x 75	Adapter	00857	1	33 places/adapter
4.0	3.0	PP Flanged Tube	03105	50	132	4,700	4,575	11 x 75	Closure	03264	50	PP Snap-On
									Adapter	00857	1	33 places/adapter
4.0	3.0	PC Flanged Tube	03104	50	132	4,700	4,575	10 x 75	Adapter	00857	1	33 places/adapter
3.0	3.0	GL Pyrex Tube	03100	50	132	4,700	4,575	10 x 75	Adapter	00857	1	33 places/adapter
1.5	1.5	High Speed Conical Microtube	314352H01	100	8	4,700	4,575	11 x 40	Adapter ⁵	00857	1	33 places/adapter
									Adapter ⁵	04871	1	33 places/adapter
—	—	Microplate (3 standard or 1 deepwell)	—	—	12	N/A	—	85 x 127; 60 mm Max height	Adapter	11779	2	Carrier for Microplates

⁵ Both adapters must be used together.

⁷ Must use with adapter cat. no. 00894. Order one for each place.

© 2010 Thermo Fisher Scientific Inc. Pyrex is a registered trademark of Corning Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

For more information, please visit
www.thermoscientific.com/rotors.

Thermo
 SCIENTIFIC