

Technical and Environmental Specifications

Below are the oven's technical and environmental specifications.

- Class 1 equipment
- Pollution degree 2
- Installation Category II
- Indoor use only in ordinary atmospheres
- Operating ambient temperatures between 15 °C and 35 °C
- Relative humidity between 15 % to 85 %
- Altitude not to exceed 2000 m operating
- G2545A is rated for mains connection to 100–120 VAC or 220–240 VAC (depending on version purchased), 50/60 Hz, 6A, with main supply voltage fluctuations up to $\pm 10\%$ of the nominal voltage

Below are the unit specifications:

Net weight	Unloaded unit weight, without rotator: 34.0 kg (75.0 lbs)
Interior dimensions	36.8 cm x 30.5 cm x 36.8 cm (14.5 in wide x 12 in deep x 14.5 in high)
Temperature	Internal operating temperature range: 5°C over ambient to 70°C. Setting resolution: 0.1°C
Rotation Speed	20 rpm