

Thermo Scientific General Purpose Rotors

TX-200 Round Buckets

TX-200 Swinging Bucket Aluminum Rotor with Round Buckets



Rotor Specifications		
Capacity (ml)	4 x 180	
Tube Dimensions (mm)	56 x 123	
Angle °	90	
Maximum Speed (rpm)	5500*	
K Factor	7920	
Net Weight (kg)	2,5	
* Maximum speed may vary by centrifuge model.		
RCF/Radius	RCF (x g)	Radius (cm)
Maximum	5.580	16,5
Minimum	2.162	6,4



TX-200 Rotor Packages

Cat. No.	Description
75003658	TX-200 Swing-Out Rotor without buckets
75003659	Round Buckets (set of 4)
75003660	Sealing Caps for Round Buckets (set of 4)
75003687	Replacement O-rings for 75003660 Sealing Caps (set of 4)

Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty.	Max. Speed (rpm)	Dims ØxL (mm)	Sealing Assembly Required:			Adapters Required:			Accessories Required:	
							Cat. No.	Qty /Pkg	Description	Cat. No.	Qty /Pkg	No. of Places		
75003800	180	180	PP Screw Cap Bottle	12	5.500	56x123	Included	12	PP Sealing	-	-	-	-	-
76009007	100	75	Glass Bottle	10	5.500	45x123	-	-	-	75003801	4	1	-	-
76009084	100	75	PP Bottle	1	5.500	45x98	-	-	-	75003801	4	1	-	-
76009095	100	75	PC Bottle	1	5.500	45 x98	-	-	-	75003801	4	1	-	-
334959	50	50	Nunc Conical	25	5.500	29.5x120	Included	25	PP Sealing	75003803	4	1	-	-
-	50	-	Universal Tube	-	-	35 x 126	-	-	-	75003802	4	1	-	-
3114-0050	50	50	Teflon Oak Ridge Tube	2	5.000	32 x 126	Included	2	PP Sealing	75003815	4	1	-	-
3139-0050	50	50	PP Oak Ridge Tube	10	5.000	32 x 126	Included	10	PP Sealing	75003815	4	1	-	-
3138-0050	50	50	PC Oak Ridge Tube	10	5.500	32 x 126	Included	10	PP Sealing	75003815	4	1	-	-
45500-30	30	25	KIMAX Glass Tube	6	-	25 x 118	-	-	-	75003805	4	2	-	-
-	25	-	Universal Tube	-	-	25,5 x 118	-	-	-	75003804	4	1	-	-
-	25	-	Round Bottom Tube	-	-	25 x 118	-	-	-	75003805	4	2	-	-
366036	15	15	Nunc Conical	50	5.500	17x120	Included	50	PP Sealing	75003771	4	5	-	-
-	15	-	Blood Collection	-	-	17 x 121	-	-	-	75003809	4	4	-	-
45500-15	15	12	KIMAX Glass Tube	50	5.500	17x102	-	-	-	75003809	4	4	-	-
-	14	-	Conical Urine Tube	-	-	17 x 114	-	-	-	75003807	4	2	-	-
3139-0010	10	8	PP Oak Ridge Tube	10	5.500	16x82	Included	10	PP Sealing	75003808	4	7	-	-
3138-0010	10	8	PC Oak Ridge Tube	10	5.500	16x82	Included	10	PP Sealing	75003808	4	7	-	-
-	10	-	Blood Collection	-	-	16x100	-	-	-	75003808	4	7	-	-
-	7	-	Blood Collection	-	-	13x100	-	-	-	75003785	4	8	-	-
-	5	-	Blood Collection	-	-	13x75	-	-	-	75003785	4	8	-	-
-	5	-	RIA Tubes	-	-	13x75	-	-	-	75003811	4	12	-	-
-	1.5/2	-	Conical Microtube	-	-	11x42	-	-	-	75003812	4	12	-	-