

**ADVANCED 1000RS ROCKING SHAKER SPECIFICATIONS**



<b>Overall dimensions (L x W x H):</b>	17 x 11 x 5" (43.2 x 27.9 x 12.7cm)
<b>Tray dimensions (L x W):</b>	12.75 x 10" (32.4 x 25.4cm)
<b>Electrical (50/60 Hz):</b>	120 volts: 25 watts 230 volts: 25 watts
<b>Fuses:</b>	5mm x 20mm, 5 amp quick acting
<b>Speed range:</b>	1 to 50rpm *
<b>Speed accuracy:</b>	±1rpm
<b>Tilt angle:</b>	0 to 15° *
<b>Timer:</b>	1 second to 9999 minutes (increased in 1 second increments)
<b>Maximum weight capacity:</b>	10lbs (4.5kg) **
<b>Controls:</b>	see page 4
<b>Tray material:</b>	aluminum
<b>Ship weight:</b>	15.5lbs (7kg)

\* Maximum speed/tilt angle may vary with heavy or unbalanced loads.  
\*\* Centered on tray.

Advanced 1000RS Rocking Shaker with glassware and plasticware