

CliniConic Fixed Angle Aluminum Rotor

Rotor Specifications

Capacity (ml)	30 x 15 ml	
Tube Dimensions (mm)	16,5 x 131	
Angle ^o	37	
Maximum Speed (rpm)	5.650	
K Factor	3954	
Net Weight (kg)	4,7	
* Max Speed may vary depending on centrifuge unit.		
RCF/Radius	RCF (x g)	Radius (cm)
Maximum	4.997	14
Minimum	3.030	8,5



Cliniconic Fixed Angle Rotor

Cat. No.	Description
75003623	CliniConic Fixed Angle Rotor

Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty.	Max. Speed (rpm)	Dims ØxL (mm)	Sealing Assembly Required:			Adapters Required:			Accessories Required:			
							Cat. No.	Qty /Pkg	Description	Cat. No.	Qty /Pkg	No. of Places	Cat. No.	Qty /Pkg	No. of Places	
-	15	15	15ml Round / Conical	-	-	16,5 x 131	-	-	-	-	-	-	-	-	-	-
-	10	-	Blood Collection	-	-	16,5 x 95	-	-	-	75003702	1	1	-	-	-	-
-	7	-	Blood Collection	-	-	16,5 x 131	-	-	-	11172595	1	1	-	-	-	-
-	5	-	Blood Collection	-	-	13 x 95	-	-	-	11172595	1	1	-	-	-	-