











Rotors	C	SC	Cat. no.	Related centrifuge	Max speed (rpm)	Max RCF (xg)
				Thermo Scientific		
Benchtop rotors						
 Fiberlite F14-6x250 LE			75003662	Sorvall™ Legend™ XT, Multifuge™ X3, SL 40F series	10,000 11,000 ^[9]	15,317 18,533 ^[9]
 Fiberlite F15-6x100y			75003698	Sorvall Legend X1/XT, Multifuge X1/X3, SL 40F series	15,000	24,652
				Sorvall ST 16/40, Megafuge™ 16/40, SL 16/40 series	13,000	18,516
 Fiberlite F13-14x50cy	■		75003661	Sorvall Legend X1, Multifuge X1 series	8,500	12,359
			75006526	Sorvall Legend XT, Multifuge X3, SL 40F series	9,250 10,000 ^[9]	14,636 17,105 ^[9]
 Fiberlite F15-8x50cy	■			75003663	Sorvall Legend XT, Multifuge X3, SL 40F series	14,500
			Sorvall Legend X1, Multifuge X1		14,000 14,500 ^[9]	22,789 24,446 ^[9]
			75006516	Sorvall Legend T, Multifuge 3 series	12,000 14,500 ^[9]	16,741 24,446 ^[9]
 Fiberlite F21-48x2			75003664	Sorvall Legend X1/XT, Sorvall ST 16/40, Multifuge X1/X3, Megafuge 16/40, SL 16/40, SL 40F series	15,200	25,055
			75006527	Sorvall Legend T, Multifuge 3 series	15,000	24,400
 Fiberlite H3-LV			75003665	Sorvall Legend X1/XT, Sorvall ST 16/40, Multifuge X1/X3, Megafuge 16/40, SL 16/40, SL 40F series	3,600	2,738

9 With 230 V centrifuge. C = Conical tubes SC = Sample containment  Biocontainment certification by Public Health England, Porton Down, UK.