



Optional, inkless, seven day, graphic, chart recorder

Control Panel:

Keyed on/off switch with set point security

Microprocessor control system with audio and visual alarms

Bright, digital temperature display

Doors:

Spring-loaded, self closing doors with 90° stay-open feature

Lockable doors for added sample security

Cabinets shown: 23.3 cu. ft. single door (RGL2304) lab refrigerator with glass door and 51.1 cu. ft. double door (RGL5004) lab refrigerator with solid doors

Thermo Scientific Revco Laboratory Refrigerators With Solid Doors

Temperature Range: +1°C to +8°C, factory preset to +4°C (auto defrost) Optional +5°C setpoint available (order option 7508TA)

Model No.	Cu. Ft. (L)	Voltage (Hz)	Certification	Plug	Cabinet	Door	Adjustable Shelves	Int. Dimensions* H x D x W inches (mm)	Ext. Dimensions* H x D x W inches (mm)	Ship Weight lb (kg)
REL404A	4.9	115 (60)	UL	NEMA 5-15	undercounter	single	3	20.0 x 20.5 x 20.0	33.4 x 26.0 x 24.0	220
REL404D	(133)	208/230 (60)	UL	NEMA 6-15				(508 x 521 x 508)	(850 x 660 x 610)	(100)
REL404V		230 (50)	CE	European						
REL404W‡		230 (50)	CE	U.K.						
Exterior height before caster installation.										
REL1204A	11.5	115 (60)	UL	NEMA 5-15	space	single	4	52.4 x 21.75 x 20.0	73.6 x 29.4 x 24.0	280
REL1204D	(326)	208/230 (60)	UL	NEMA 6-15	saver			(1331 x 553 x 508)	(1869 x 748 x 610)	(127)
REL1204V		230 (50)	CE	European						
REL1204W‡		230 (50)	CE	U.K.						
REL2304A	23.3	115 (60)	UL	NEMA 5-15	upright	single	4	58.0 x 29.0 x 24.0	79.2 x 37.2 x 28.0	370
REL2304D	(659)	208/230 (60)	UL	NEMA 6-15				(1473 x 737 x 610)	(2012 x 945 x 711)	(168)
REL2304V		230 (50)	CE	European						
REL2304W‡		230 (50)	CE	U.K.						
REL3004A	29.2	115 (60)	UL	NEMA 5-15	upright	single	4	58.0 x 29.0 x 30.0	79.2 x 37.2 x 34.0	435
REL3004D	(826)	208/230 (60)	UL	NEMA 6-15				(1473 x 737 x 762)	(2012 x 945 x 864)	(197)
REL3004V		230 (50)	CE	European						
REL3004W‡		230 (50)	CE	U.K.						
REL5004A	51.1	115 (60)	UL	NEMA 5-15	upright	double	8	58.0 x 29.0 x 52.5	79.2 x 37.2 x 56.5	585
REL5004D	(1447)	208/230 (60)	UL	NEMA 6-15				(1473 x 737 x 1334)	(2012 x 945 x 1435)	(265)
REL5004V		230 (50)	CE	European						
REL5004W‡		230 (50)	CE	U.K.						
REL7504A	78.8	115 (60)	UL	NEMA 5-15	upright	triple	12	58.0 x 29.0 x 81.0	79.2 x 37.2 x 85.0	805
REL7504D	(2231)	208/230 (60)	UL	NEMA 6-15				(1473 x 737 x 2057)	(2012 x 945 x 2159)	(365)
REL7504V		230 (50)	CE	European						
REL7504W		230 (50)	CE	U.K.						

► **Would you like to customize your solid door refrigerator? See page 17 for options and accessories.**

* Published interior volumes are nominal; useable cu. ft. (L) capacities may be offset by shelves, drawers, blower coils, or cabinet protrusions. For specific cabinet dimensions and allowances, or for verification of clearance tolerances, contact our sales representative. Product performance may be affected by ambient temperature conditions, humidity conditions, sample loading preferences, and frequency of door openings. Contact our sales representative for specific inquiries regarding your application. 4.7 cu. ft. exterior dimensions are without casters.

‡ Due to differences in building wiring regulations in Australia and New Zealand, the current draw for the W-voltage units rated less than 10A may utilize a 10A circuit breaker. This rating only applies to Australia and New Zealand. Confirm the amperage rating on the data plate. Contact your local sales representative for more information.