

5.0 SPECIFICATIONS

	SPD1010	SPD2010
Capacity:	Moderate	Large
Refrigerated trap:	-50 °C (approx) 4 liter CFC-free	-50 °C (approx) 4 liter CFC-free
Vacuum Pump:		
Displacement (50/60 Hz)	30/36 liters/min	30/36 liters/min
Maximum vacuum	7 Torr	7 Torr
Vacuum Level Control:	30 to 5.1 torr in 0.1 increments	
Vacuum Ramp Setting:	1 to 5 .	
Size (W X D X H):	25 in x 26 in x 15 in 62 cm x 65 cm x 37 cm	31 in x 27 in x 19 in 77 cm x 67 cm x 47 cm
Weight:	152 lbs 69 kg	182 lbs 83 kg
Operative Power:	115VAC; 60 Hz; 12 amps 230 VAC; 50 Hz; 6 amps	220 VAC; 60 Hz; 8 amps 230 VAC; 50 Hz 8 amps
*Dependent upon ambient temperature, line voltage fluctuation, and load capacity.		