


MANUAL ANALOG BENCHTOP AUTOCLAVES

Designed to autoclave glass, wrapped, and unwrapped tools.
Not suitable for liquids (see liquid autoclaves on pages 8 to 12).

- 60 minute cycle timer
- Temperature selection from 100° to 134°C (212° to 273°F)
- Combination pressure and temperature gauge
- Includes trays



	1730M	2340M	2540M	3870M
Chamber Volume (L)	7.5	19	23	85
Max Available Height (in)	5.25"	6.75"	8.3"	12.1"
Chamber Dimension (ØxD)	7" x 12"	9" x 18"	10" x 18"	15" x 27"
External Dimension (WxHxD)	17" x 12.5" x 16"	20" x 15" x 20"	20" x 15" x 20"	26" x 21" x 34"
Power	120V	120V	120V	230V, 1 phase (60 Hz) NEMA 6-15 
Utilities	No Plumbing Required			
P/N	023210002	023210100	023210304	023210550
Autoclave Stand (PN 023212609)				Free with Autoclave
Table with SS top (30" x 36") & Castors	023212601	023212601	023212601	023212601

2 YEAR WARRANTY

ACCESSORIES

Printer Paper		023212200
Chamber Brite		023212805
Mobile Lifting Device for 5075ELVC & 3870ELV-WR		023250016
Biohazard Air Filtration System - Factory Installed		023250011
5075ELVC Large SS Wire Basket	023211790	Model 3850 SS Tray Large with no holes 023213583
5075ELVC Small SS Wire Basket	023211795	Model 3850 SS Tray Small with holes 023213591
2840ELV -High SS Wire Basket	023211805	Model 3870 SS Tray Small with no holes 023213672
2840ELV -Lower SS Wire Basket	023211806	Model 3850 SS Tray Small with no holes 023213753
2840ELV -Lower SS Solid Basket	023211807	Model 2840ELP Solid Tray 023290011
3850/70ELV -Large SS Wire Basket	023211840	Pouch Rack for 2340, 2540, 3545 023212706
3850/70ELV -Small SS Wire Basket	023211850	Pouch Rack for 3870 023212805
3850/70ELV -Small SS Solid Basket	023211860	Model 2840ELP Wire SS Basket, Wire 400x180x160mm 023211809
3850/70ELV -Large SS Solid Basket	023211870	Model 5075ELP Wire SS Basket 715x380x250mm 023211891
3850/70ELP SS Wire Basket	023211871	Model 38xx Wire SS Basket 500x260x225mm 023211892
5075ELVC Large SS Solid Basket	023211890	Model 2840ELP SS Container 400x190x160mm 023211808
5075ELVC Small SS Solid Basket	023211895	