

UNIT SPECIFICATIONS

These units are 115 volt or 220 volt. Please refer to the unit data plate for its individual specification.

Weight	Shipping	Net
2300	280 lbs.	221 lbs.
2310	280 lbs.	257 lbs.
2350	600 lbs.	442 lbs.
2400	280 lbs.	221 lbs.
2450	800 lbs.	442 lbs.

Dimensions	Exterior WxDxH (IN.)	Interior WxDxH (IN.)
2300	26.5x26x40	20.5x20x28
2310	26.5x25.25x39.75	20.5x20x25.75
2350	26.5x26x80	20.5x20x28
2400	26.5x26x40	19X20X25.5
2450	26.5x26x80	19X20X25.5 (x2)

Capacity	Cubic Feet
2310	2
2300	6.7
2350	6.7 (x2)
2400	6.7
2450	6.7 (x2)

Temperature	Range	Uniformity	CO2
2300	Amb. 8°C to 60°C	$\pm 0.1^{\circ}\text{C}$ @ 37°C	0 – 20%
2310	Amb. 8°C to 60°C	$\pm 0.1^{\circ}\text{C}$ @ 37°C	0 – 20%
2350	Amb. 8°C to 60°C	$\pm 0.1^{\circ}\text{C}$ @ 37°C	0 – 20%
2400	Amb. 8°C to 60°C	$\pm 0.1^{\circ}\text{C}$ @ 37°C	0 – 20%
2450	Amb. 8°C to 60°C	$\pm 0.1^{\circ}\text{C}$ @ 37°C	0 – 20%

Model	Amp	Voltage	Hertz
2300	5 amp	120V	60 Hz

Model	Amp	Voltage	Hertz
2300-2	3.5 amp	220V	50/60 Hz
2310	5 amp	120V	60 Hz
2310-2	3.5 amp	220V	50/60 Hz
2350	5 amp	120V	60 Hz
2350-2	3.5 amp	220V	50/60 Hz
2400	5 amp	120V	60 Hz
2400-2	3.5 amp	220V	50/60 Hz
2450	5 amp	120V	60 Hz
2450-2	3.5 amp	220V	50/60 Hz