

## SPECIFICATIONS - HIGH SPEED MICROPLATE SHAKER



<b>Overall dimensions (L x W x H):</b>	15.5 x 12 x 12.8" (39.4 x 30.5 x 32.5cm)
<b>Tray dimensions (L x W):</b>	12.25 x 11" (31.1 x 27.9cm)
<b>Electrical (50/60 Hz):</b>	120 V - 1 amps, 120 watts 230 V - 1 amp, 120 watts
<b>Fuses:</b>	5mm x 20mm, 5 amp quick acting (2 required)
<b>Speed range:</b>	600 to 2500rpm
<b>Timer:</b>	1 second to 160 hrs
<b>Orbit:</b>	.014" (3.6mm)
<b>Controls:</b>	See page 5
<b>Ship weight:</b>	50lbs (22.7kg)
<b>Capacity:</b>	7 lbs (3.2kg) 48 microplates in stacks up to 5" high