

## VWR<sup>®</sup> Mini Vortex Mixer



### VWR<sup>®</sup> Mini Vortex Mixer

The VWR Mini Vortex Mixer is the ideal vortex mixer for any laboratory. Its simple "Touch & Go" design and function make it easy to use. Small and compact to conserve precious bench space, yet powerful enough to deliver robust vortex mixing at 3,000 rpm.

- Compact and sleek design
- 3,000 rpm for powerful vortex mixing
- Touch activated for simple use
- 2 year warranty

Capacity	Electrical	Dimensions	Speed	Orbit	Maximum Power	Shipping Weight	Cat. No.
80 mL	110V, 60 Hz	96 x 85 mm (3 <sup>13</sup> / <sub>16</sub> x 3 <sup>3</sup> / <sub>8</sub> "	3,000 rpm	4.5 mm ( <sup>3</sup> / <sub>16</sub> "	5 W	0.55 Kg (1.2 lbs.)	10153-688

**VWR**  <sup>®</sup> **1.800.932.5000**  
We Enable Science **vwr.com**

Prices and product details are current when published; subject to change without notice. [ Certain products may be limited by federal, state, provincial, or local regulations. ] VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. [ VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). ] Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners and additional disclaimers. ©2015 VWR International, LLC. All rights reserved.