

Tanks – 1295PC	7930510	7930510
----------------	---------	---------

UNIT SPECIFICATIONS

Weight	Shipping	Net
1225, PC	21 lbs.	14 lbs
1235, PC	26.5 lbs.	18 lbs.
1245, PC	29 lbs.	25 lbs.
1255, PC	29 lbs.	32 lbs.
1295PC	34.3 lbs.	28.1 lbs.

Dimensions	Exterior WxDxH (in.)	Interior WxDxH (in.)
1225, PC	14 x 12.5 x 9.8	11.8 x 5.8 x 6
1235, PC	14 x 19.3 x 9.5	11.8 x 12.8 x 5.8
1245, PC	14 x 26 x 10	11.5 x 19.5 x 6
1255, PC	23.5 x 19.3 x 9.3	26.5 20 x 20.5
1295PC	16 x 19 x 16	12 x 14 x 11

Capacity	Liters
1225, PC	6
1235, PC	14
1245, PC	20
1255, PC	6 & 14
1295PC	30

Temperature	Range	Uniformity
1225, PC	Amb. 5° to 100°C	±0.2°C @ 37°
1235, PC	Amb. 5° to 100°C	±0.2°C @ 37°
1245, PC	Amb. 5° to 100°C	±0.2°C @ 37°
1255, PC	Amb. 5° to 100°C	±0.2°C @ 37°
1295PC	Amb. 5° to 100°C	±0.3°C @ 37°